



Lecture 9
27/01/2023
By Prof Manish Kumar

Observing Microorganisms



Units of Measurement

- $1\ \mu\text{m} = 10^{-6}\ \text{m} = 10^{-3}\ \text{mm}$
- $1\ \text{nm} = 10^{-9}\ \text{m} = 10^{-6}\ \text{mm}$
- $1000\ \text{nm} = 1\ \mu\text{m}$
- $0.001\ \mu\text{m} = 1\ \text{nm}$




E. coli
bacteria

SEM

1 μm

Microscopy: The Instruments

- A simple microscope has only one lens




(b) Microscope replica

Light Microscopy


- Use of any kind of microscope that uses visible light to observe specimens
- Types of **light microscopy**
 - Compound light microscopy
 - Darkfield microscopy
 - Phase-contrast microscopy
 - Differential interference contrast microscopy
 - Fluorescence microscopy
 - Confocal microscopy

The Compound Light Microscope



Compound Light Microscopy

- In a **compound microscope**, the image from the objective lens is magnified again by the ocular lens
- **Total magnification** = objective lens \times ocular lens



(b) The path of light (bottom to top)

Figure 3.1b


Compound Light Microscopy

- **Resolution** is the ability of the lenses to distinguish two points
- A microscope with a resolving power of 0.4 nm can distinguish between two points ≥ 0.4 nm
- Shorter wavelengths of light provide greater resolution

Compound Light Microscopy



- The **refractive index** is a measure of the light-bending ability of a medium
- The light may bend in air so much that it misses the small high-magnification lens
- Immersion oil is used to keep light from bending

Refraction in the Compound Microscope



Brightfield Illumination



- Dark objects are visible against a bright background
- Light reflected off the specimen does not enter the objective lens



LM 20 μm

Darkfield Illumination



- Light objects are visible against a dark background
- Light reflected off the specimen enters the objective lens



LM 20 μm

Phase-Contrast Microscopy

- Accentuates (emphasize) diffraction of the light that passes through a specimen
- The wave nature of light rays.
- Light rays can be in phase (their crests and trough match) or out of phase
- Light rays (direct and diffracted/refracted) interact to produce reinforcement (relative brightness) or interference (relative darkness)



LM 20 μm


Differential Interference Contrast Microscopy

- Accentuates diffraction of the light that passes through a specimen; uses two beams of light
- Employs polarizer in the condenser.
- The polarized light is passed through prism and generates two beams
- Image is brightly colored and appears 3-D



Fluorescence Microscopy

- Uses UV light
- Fluorescent substances absorb UV light and emit visible light
- Cells may be stained with fluorescent dyes (fluorochromes)



(b)

LM

5 μ m

Lecture 10
30/01/2023
By Prof Manish Kumar


Confocal Microscopy

- Cells stained with fluorochrome dyes
- Short wavelength (blue) light used to excite the dyes
- The light illuminates each plane in a specimen to produce a three-dimensional image
 - Up to 100 μm deep



Scanning Acoustic Microscopy (SAM)

- Measures sound waves that are reflected back from an object
- Used to study cells attached to a surface (cancer cells, artery plaque)
- Resolution $1\text{ }\mu\text{m}$




SAM | 170 μm

Electron Microscopy

- Uses electrons instead of light
- The shorter wavelength of electrons gives greater resolution

Transmission Electron Microscopy (TEM)

- Ultra thin sections of specimens
- Specimens are placed on copper mesh grids
- Beam of electrons passes through specimen, then an electromagnetic lens, to a screen or film
- Specimens may be stained with heavy metal salts
- Magnify 10000 to 100000X.



Transmission Electron Microscopy (TEM)


- 10,000–100,000 \times ; resolution 2.5 nm



TEM 15 μ m

Scanning Electron Microscopy (SEM)

- Overcomes problem of sectioning specimen
- An electron gun produces a beam of electrons that scans the surface of a whole specimen
- Secondary electrons emitted from the specimen produce the image



Scanning Electron Microscopy (SEM)


- 1,000–10,000 \times ; resolution 20 nm
- Useful for surface structures of intact cells and viruses.



SEM 15 μ m

Scanned-Probe Microscopy

- **Scanning tunneling microscopy (STM)** uses a metal probe to scan a specimen
- Resolution 1/100 of an atom (thus greater than EM)



(a)


RecA protein of *E. coli*



Lecture 11
31/01/2023
By Prof Manish Kumar

Scanned-Probe Microscopy

- **Atomic force microscopy (AFM)** uses a metal- and-diamond probe inserted (stylus) into the specimen.
- Produces three-dimensional images.



(b)

Perfringoglysin O toxin





Figure 3.11b

Preparing Smears for Staining

- **Staining:** Coloring the microbe with a dye that emphasizes certain structures
- **Smear:** A thin film of a solution of microbes on a slide
- A smear is usually **fixed** to attach the microbes to the slide and to kill the microbes


Preparing Smears for Staining

- Live or unstained cells have little contrast with the surrounding medium. Researchers do make discoveries about cell behavior by observing live specimens.



5 μm

LM



10 μm

SEM

Preparing Smears for Staining

- Stains are salts composed of a positive and negative ion
- In a **basic dye** (methylene blue, malachite green and safranin), the chromophore is a cation
- Bacteria are slightly negatively charged at pH7.
- In an **acidic dye**, the chromophore is an anion
- Acidic dye stain the background instead of the cell is called **negative staining**

Simple Stains

- **Simple stain:** Use of a single basic dye. Eg. Methylene blue, carbolfuchsin, crystal violet, and safranin.
- A **mordant** may be used to hold the stain or coat the specimen to thicken it for easier visibility

Differential Stains

- Used to distinguish between bacteria
 - Gram stain
 - Acid-fast stain


Gram Stain

- Classifies bacteria into gram-positive or gram-negative
 - Gram-positive bacteria tend to be killed by penicillin and detergents
 - Gram-negative bacteria are more resistant to antibiotics

Gram Stain

	Color of Gram-positive cells	Color of Gram-negative cells
Primary stain: Crystal violet	Purple	Purple
Mordant: Iodine	Purple	Purple
Decolorizing agent: Alcohol-acetone	Purple	Colorless
Counterstain: Safranin	Purple	Red

Micrograph of Gram-Stained Bacteria



(b)

LM 5 μm

Acid-Fast Stain


- Stained waxy cell wall is not decolorized by acid-alcohol
- *Mycobacterium*
- *Nocardia*

Acid-Fast Stain

	Color of Acid-fast	Color of Non-Acid-fast
Primary stain: Carbolfuchsin	Red	Red
Decolorizing agent: Acid-alcohol	Red	Colorless
Counterstain: Methylene blue	Red	Blue

Acid-Fast Bacteria

M. leprae




LM 5 μ m

Special Stains

- Used to distinguish parts of cells
 - Capsule stain
 - Endospore stain
 - Flagella stain

Negative Staining for Capsules

- Cells stained
- Negative stain




(a) Negative staining

LM 5 μ m

Endospore Staining

- Primary stain: Malachite green, usually with heat
- Decolorize cells: Water
- Counterstain: Safranin




(b) Endospore staining

LM 5 μ m

Flagella Staining

- Mordant on flagella
- Carbofuchsin simple stain



(c) Flagella staining

LM 5 μ m