


5. Calculate the distance between the donor and acceptor molecules from the given Förster resonance energy transfer data if the Förster distance for the given FRET pair is 65 \AA . {4 marks}



6. You are studying the fluorescence quenching of a peptide containing a single Trp. The acrylamide quenching data recorded at emission maximum ($\lambda_{\text{ex}} = 295 \text{ nm}$) is shown below.

Concentration of acrylamide (M)	Fluorescence intensity (AU)	Lifetime (ns)
0	1000	17.6
0.05	216	3.26
0.1	94	1.8
0.2	43	0.95
0.3	27	0.64
0.4	15	0.49
0.5	7	0.39

- Trp absorption maximum is near 280 nm , what could be the reason of using 295 nm excitation wavelength? {1 mark}
- Draw a neat Stern-Volmer plot. {3 marks}
- Is the quenching static or dynamic or both? {1 mark}
- Determine the quenching constant(s) (K_0 , K_s , or both: depends on your answer in part 'b') {3 marks}

APPENDIX

$$E = \frac{K_0}{K_0 + I^0}$$

$$\frac{F_0}{F} = 1 + (K_0 + K_s) [Q] + K_0 K_s [Q]^2$$

$$\frac{F_0}{F} = 1 + \tau F_0 k_q [Q] = 1 + 6D\tau$$

$$r = \frac{I_0 - I_\infty}{I_\infty + 2I_\infty}$$

$$I(t) = I_0 e^{-t/\tau}$$

INDIAN INSTITUTE OF TECHNOLOGY GUWAHATI

MID-SEMESTER EXAMINATION, 2023

BT 624: Fluorescence Techniques in Biotechnology

Date: March 03, 2023

Time: 2 - 4 P.M.

Maximum marks: 30

Name: Bhargya Guddake

Roll No. 200166020

Instructions

1. Write your name and Roll No. on the answer sheet.
2. The question paper carries 6 questions that span 2 pages.
3. Some formulas are given in the Appendix. You can use them if required, otherwise ignore them.
4. You will be given a Graph sheet. If you need more than one, please ask.

Attempt all questions

1. a. Draw a labeled Jablonski diagram showing absorption, internal conversion, fluorescence, and phosphorescence. {3 marks}
b. Arrange the processes mentioned in part (a) with decreasing rates (increasing lifetime). {1 mark}
2. Explain briefly the following terms: {6 marks}
 - a. Fluorescence
 - b. Phosphorescence
 - c. Singlet state
 - d. Triplet state
 - e. Quantum yield
 - f. Intersystem crossing
3. Using an appropriate diagram, explain the Franck-Condon principle and mirror image rule for fluorescence emission. {4 marks}
4. The intrinsic fluorescence anisotropy, r_0 of a fluorophore is 0.36. Calculate the hydrodynamic volume of the fluorophore from the given data if the observed anisotropy in water at 20 °C is 0.1.
Coefficient of viscosity of water at 20 °C = 0.1 Pa.s
Fluorescence lifetime: 6 ns
Rotational correlation time = $\frac{\eta V}{RT}$,
where, η is the coefficient of viscosity, V is the hydrodynamic volume of the rotating unit, R is the Universal gas constant ($8.314 \text{ J mol}^{-1} \text{ K}^{-1}$ or $8.314 \text{ Pa.m}^3 \text{ mol}^{-1} \text{ K}^{-1}$) and T is the temperature. {4 marks}